

Want to know more?
Give Lennox a call on 01 455 2201

memmert
Experts in Thermostatics

CO₂ Incubators **INCO 2**



*Safe and precise control of CO₂,
humidity and temperature*

*Humidity control from 88%-97% rh – optional active
humidifying system adjustable from 40% to 97% rh*

*Working chamber including ventilation system, water trays as
well as all sensors sterilisable without problems*

Triple precision

For safe processes

Contamination? Not with us!

>>>> www.memmert.com

Available in Ireland from



Want to know more?
Give Lennox a call on 01 455 2201



Stainless steel – the noble material

The new INCO 2 – this is fine-tuning to perfection: for absolutely reliable experimental procedures with a maximum of security and operating convenience.

The inner chamber of the INCO 2 consists exclusively of high-grade and fully recyclable stainless steel. The additional electropolishing of the chamber ensures particularly smooth and hygienic surfaces and therefore the greatest possible reduction in contamination risk. The outer casing is made mainly from sturdy textured stainless steel (except rear of zinc-plated steel) and is in accord with Memmert's holistic quality philosophy.

Bacteria-free through sterilisation function – desinfection is not sufficient!

Hygiene is fundamental in high-sensitivity applications with organic materials. To avoid cross-contamination the entire chamber, including the ventilation system, the water tray and all sensors can be sterilised at 160°C in a 4-hour programme. The infrared CO₂ sensor was specially designed by Memmert engineers to withstand without damage the 160°C sterilisation temperature over the entire 4-hour period. For safety reasons this function can only be operated by means of the Memmert STERICard in order to prevent unintentional activation.

The professional CO₂ – working for man and nature

CO₂-atmosphere up to 20%, accurate temperature control up to 45°C and humidity limit control up to 97% rh without condensation combine to create inside the INCO 2 a controlled and physiologically ideal environment for cell and tissue cultures growing in vitro.

3 sizes:
108 - 153 - 246 litre chamber volume



Want to know more?
Give Lennox a call on 01 455 2201

Uniformity inside the chamber

Chamber heating from all six sides, together with electronic humidity limit control, is mainly responsible for preventing condensation in the chamber. An aluminium heat-conducting jacket helps to ensure optimal temperature distribution and in addition serves as heat store in case of a temporary supply outage.

Despite the appreciably higher density of CO₂ compared with air, the turbulence-free chamber ventilation provides a uniformly homogeneous atmosphere.



A Memmert special: the chamber is made to 100% of high-grade stainless steel mat. 1.4301 (ASTM 304) and is additionally electropolished for maximum protection against dangerous contamination.


A Memmert special: the all-round heating is located below deep-drawn ribs which carry the anti-tipping stainless steel shelves and at the same time ensure optimal and particularly gentle heat transfer to the incubator load.

Active humidity control

In addition to CO₂ content, a stable and accurately adjustable humidity plays a decisive role in creating physiologically ideal ambient conditions for the incubator load.

Memmert offers the INCO 2 in two variants. As standard it is equipped with a humidity limit control by which the maximum relative humidity created by water dishes in the chamber can be reduced from 97% down to 88%. At extra charge an active humidification system can be incorporated, with controlled evaporation of water through a steam generator and sterile dry steam mixed into the air flow. This system offers an extended control range of 40% to 97% rH.

Your advantage: more useful volume, reduced contamination risk and short recovery time after opening the door.



The difference in humidity recovery time between humidity limit control and active humidity control (example: INCO 2/108).



Want to know more?
Give Lennox a call on 01 455 2201



There is no such thing as a little reliability!

The new INCO 2 is a model of reliability and accuracy. Innovations in control engineering ensure maximum security and reliability. And surprise: this breadth of functionality is obtained with even more convenience and ease of operation!

Convenience: made for you

The INCO 2 is easy to use:

- Clear, easy-to-clean underglass function display
- Memmert's unique feature: the push-turn control (patented) for intuitive selection from the entire menu
- Feet with height adjustment as standard
- Auto start-up programme
- Automatic gas cylinder changeover
- Quick-release connectors on CO₂ intake, also for water supply on active humidity control
- Fully insulated stainless steel door and additional glass door

Accuracy: for regulated procedures


Technical features for faultless processes:

- Internal ventilation for turbulence-free and uniform distribution of CO₂, humidity and temperature
- Multifunctional fuzzy-supported control for exact setting and maintenance of temperature, humidity and CO₂ content
- Calibration facility for temperature, CO₂ and humidity on the actual controller
- Option (extra charge): up to three additional and freely positionable Pt100 sensors for indication on the display and recording in the document store
- Digitised NDIR CO₂ sensor with cyclic auto-zero function
- Compensation for atmospheric pressure errors through internal barometric pressure measurement

Documentation: for controlled quality

The basic outfit for professional quality assurance:

- Standard "Celsius" software for programming and documentation, also optional FDA-conforming software (extra charge)
- Internal ring memory for uninterrupted, long-term documentation fully protected against manipulation (approx. 3 months).
- Parallel printer interface for printouts of thermostating processes (USB printer via converter possible)
- USB interface for programming and storage (Ethernet at extra charge)




Security: of course!

Still more functions with no error:

- Integral auto-diagnostic system with visual and audible error indication
- Triple temperature monitor with mechanical temperature limiter (TB), electronic adjustable temperature monitor (TWW), as well as the Memmert ASF (Auto-Safety Function)
- Two high-grade platinum sensors communicate with each other like a climbing team on a rope and secure uninterrupted faultless temperature control. The high-grade 4-wire circuit ensures precise transmission of the measurements
- Audible signal on over- and undertemperature, over- and under-CO₂, open door and empty gas cylinder
- Option (extra charge): volt-free contacts on "combination error" and "setpoint reached"
- Option (extra charge): personal User-ID-Card for protection against unauthorised operation

Want to know more?
Give Lennox a call on 01 455 2201



Timer module


- 1 Time indication (real-time)
- 2 Text messages

Temperature module

- 11 Heating
- 12 Sterilisation operation
- 13 Setpoint/actual temperature
- Temperature adjustment:
20 – 50°C (control range from 8°C above ambient temperature up to 50°C)
- Variation (time):
±0.1°C max.
- Uniformity (spatial):
±0.3°C max.

CO₂ module

- 24 Gas cylinder 1 (active)
- 25 Gas cylinder 2
- 26 CO₂ setpoint/actual
- CO₂-concentration selection 0-20% (adjustable in 0.1% steps)
- Variation (time)
±0.1% CO₂ max.
- Uniformity (spatial)
±0.3% CO₂ max.



Operating mode

- 3 Normal operation (active)
- 4 Weekly timer*
- 5 Printer
- 6 Configuration
- 7 Data manipulation prevented through optional User-ID-Card (extra charge)
- 8 SET key
- 9 Push-turn control
- 10 Chip card reader for STERICard and optional User-ID-Card

Monitor module

- 14 Visual alarm
- 15 Alarm limit (heating switch-off temperature)
- 16 Sounder
- 17 Low alarm limit
- 18 Automatic alarm limit (ASF)
- 19 Temperature limiter
- 20 High alarm limit

Audible and visual alarm on temperature and CO₂ out of limit, on door open longer than 2 minutes, and on other errors

Humidity module

- 21 Setpoint/actual humidity
- 22 Water level in storage tank
- 23 Steam generation
- Humidity selection 40-97% rh with active humidity control
- Humidity selection 88-97% rh with humidity limit control (selected in 1% steps)
- Variation (time)
±1% rh max.

* Weekly timer programmable with one ON and OFF time per weekday; additional group function

Want to know more?
Give Lennox a call on 01 455 2201

Technical data, models and accessories for CO₂ Incubators INCO 2


Interfaces:



optional
(extra charge): Ethernet LAN



Dimensions of CO₂ Incubators
(see table below)



Model sizes		INCO 2	108	153	246
Stainless steel interior (deep-drawn, electropolished)	Volume	approx. l	108	153	246
	Width / Height / Depth (less 25 mm for air duct in the middle of the back wall) (A) / (B) / (C)	mm	560 / 480 / 400	480 / 640 / 500	640 / 640 / 600
	Provision for sliding stainless steel shelves or wire grid shelves half width / full width	number	- / 4	- / 6	2 x 6 / 6
Stainless steel exterior (rear zinc-plated steel)	Width (D)	mm	710	630	790
	Height (variable through adjustable feet) (E)	mm	778	938	938
	Depth (without door handle, depth of handle 38 mm) (F)	mm	550	650	750
	Fully insulated, heated stainless steel door		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Extra internal glass door		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Ventilation	Uniform atmosphere and temperature distribution through enclosed non-turbulent ventilation system, fully covered by the sterilisation process		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Temperature	Electronic microprocessor temperature controller with Pt100 and auto-diagnostic system		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Temperature sensors Pt100 Class A in 4-wire circuit for uninterrupted operation on failure of one Pt100 with warning indication		double	double	double
	Temperature range (during sterilisation the temperature is fixed at 160 °C – set value)	° C	from 20 ¹⁾ up to 50	from 20 ¹⁾ up to 50	from 20 ¹⁾ up to 50
	Temperature fluctuations with time (to DIN 12 880: 2007-05)	° C	≤ ± 0,1	≤ ± 0,1	≤ ± 0,1
	Temperature variation in chamber at 37 °C (to DIN 12 880: 2007-05)	° C	≤ ± 0,3	≤ ± 0,3	≤ ± 0,3
Sterilisation	STERICard for automatic chamber sterilisation cycle 4 h at 160 °C (not for sterilising the load!)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
CO ₂	Digital electronic CO ₂ control with autozero, NDIR system and integrated gas cylinder change-over, with auto-diagnostic system and acoustic fault indication		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Setting accuracy	CO ₂	0,1 %	0,1 %	0,1 %
	Adjustment range	CO ₂	0 up to 20 %	0 up to 20 %	0 up to 20 %
Humidity	Capacitive humidity sensor (sterilisable)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Standard water dishes	number	1	1	2
	Adjustable humidity limit control (88 – 97%) incl. digital indication and auto-diagnostic system with visual and acoustic fault indication (air supply via sterile filter) ensures rapid reaching of set humidity and short recovery times while avoiding condensate formation		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Monitor	Microprocessor temperature monitor acting as overtemperature protection, with Pt100 incorporating fault diagnostics with visual and audible alarm		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Digital over- and undertemperature monitor		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Temperature monitoring band automatically linked to the setpoint (ASF)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Relay for reliable heating cut-off in case of fault		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Mechanical temperature limiter (TB)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Audible alarm: Over- and undertemperature, over- and under-CO ₂ , open door and empty gas cylinder		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Timer functions	Real-time/weekly programmer with group function (e.g. Monday-Friday)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Documentation	Internal log memory 1024 kB as ring memory for all setpoints, actual values, errors, settings with real-time and date; capacity approx. 3 months at 1 min intervals		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Parallel printer interface for printing logging files, suitable for all PCL3-compatible ink jet printers (USB available via converter, see accessories)		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	"Celsius" ²⁾ software for control and documentation of temperature, CO ₂ and relative humidity		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Setup	Calibration (no separate PC required): Temperature: 3-point calibration on controller. CO ₂ : 3-point calibration at 5%, 10% and 15%, Auto-zero-function of NDIR CO ₂ -sensor upon switching on and cyclically every 24 h. Humidity: 2-point calibration at 20% and 90%		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Setting of language for dialogue and display D / UK / E / F / I		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Want to know more?

Give Lennox a call on 01 455 2201

Model sizes		INCO 2	108	153	246
Further data	Electrical supply ³⁾	V / Hz	230 / 50/60	230 / 50/60	230 / 50/60
	Electrical load	approx. W	1000	1500	2000
	Net weight / Gross weight	approx. kg	70 / 78	82 / 114	110 / 160
Standard accessories	2 connecting hoses with two quick-release couplings to 2 gas cylinders (automatic gas cylinder change over)		□	□	□
	Stainless steel sliding shelves (full width)	number	2 □	3 □	3 □
	Shelf or wire grid shelf width (half width / full width)	approx. mm	- / 556	- / 480	300 / 640
	Shelf depth	approx. mm	361	400	500
	Stainless steel water dishes, electropolished	number	1 (full width)	1 (full width)	2 (half width)
	Width of dish	approx. mm	535 (full width)	455 (full width)	280 (half width)
	Depth of dish / height of dish	approx. mm	250 / 40	350 / 40	450 / 40
	Works calibration certificate (test point chamber centre at 37 °C)		□	□	□
Standard version	Incubator INCO 2		INC 108	INC 153	INC 246
Options	4-part partitioning of interior with 4-part gas baffle (replacement of 3 full-width shelves by 6 half-width wire grid shelves)		–	–	K4
	Active microprocessor control for humidifying and dehumidifying (40 – 97 % rh), incl. digital indication and auto-diagnostic system ensures even more rapid reaching of set humidity and very short recovery times while avoiding condensate formation; humidity supply with distilled water (from an external tank) by a self-priming pump; integral bacteria block by generating hotsteam, dehumidifying via sterile filter; (standard humidity limit control and water dishes are omitted). Can only be ordered together with new incubator (no retrofitting possible)		K7	K7	K7
	Additional Pt100 temperature sensor, positioned flexibly in chamber or load, for local temperature measurement (up to 3 additional sensors are possible). The measured temperature can, if required, be indicated on the multifunction display, recorded in the integral ring store, and can be documented via the "Celsius" ²⁾ software or on an attached printer		H8(x)	H8(x)	H8(x)
	Extra stainless steel shelf, half width		–	–	E6(x)
	Extra stainless steel shelf, full width		E7(x)	E7(x)	E7(x)
	Extra wire grid shelf, full width		E3(x)	E3(x)	E3(x)
	Extra water dish		E2(x) (full width)	E2(x) (full width)	E2(x) (half width)
	Stacking version for 2 units of equal size (bottom unit modification)		G3	G3	G3
	Subframe (622 mm high) ⁴⁾		G5	G5	G5
	Subframe (130 mm high for 2 stacked incubators) ⁴⁾		G7	G7	G7
	STERICard (additional or as replacement) for automatic chamber sterilisation cycle (not for sterilising the load!)		V9	V9	V9
	Oven-linked authorisation card (User-ID-Card) prevents undesired manipulation by unauthorised third parties		V1	V1	V1
	Computer interface RS485 (for networking a maximum of 16 ovens) instead of USB		V2	V2	V2
	Interface RS232 instead of USB		W6	W6	W6
	Interface Ethernet instead of USB incl. software "Celsius Ethernet-Edition"		W4	W4	W4
	Parallel/USB converter cable with integrated power supply unit to connect HP printers with USB interface to MEMMERT units. Compatible with USB 1.1 and USB 2.0		W1	W1	W1
	Documentation package consisting of parallel USB converter cable including PLC3-compatible HP colour inkjet printer with USB interface (HP Deskjet 6940 or successor) for direct connection of printer to Memmert unit		W2	W2	W2
	USB connection cable for computer interface		W7	W7	W7
	Connection cable for computer interface RS232 according to DIN 12 900-1		V6	V6	V6
	Potential-free contact (24 V/2 A) with socket to NAMUR NE 28 for external monitoring (setpoint of temperature and CO ₂ is reached)		H5	H5	H5
	Potential-free contact (24 V/2 A) with socket to NAMUR NE 28 for combination error message of temperature and CO ₂ controller (e. g. supply failure, sensor fault or fuse)		H6	H6	H6
	Pressure reducing valve to DIN 8546, incl. gas cylinder monitor		H7	H7	H7
	Works calibration certificate for 5%, 10% and 15% CO ₂ (measured at 37 °C)		Z5	Z5	Z5
IQ check list with works test data for incubator as support for validation by customer		Q1	Q1	Q1	
OQ check list including one ⁵⁾ free-selectable temperature distribution survey for 27 measuring points to DIN 12880: 2007-05 with works test data for 5% CO ₂ and 90% humidity as support for validation by customer		Q2	Q2	Q2	
Software "Celsius FDA Edition" ⁶⁾ for up to 16 units. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA)		Q3	Q3	Q3	
Integration of additional units (up to max. 16 units) into an already existent FDA-software licence		Q4	Q4	Q4	
Start-up of INCO incubators and brief training ⁷⁾ (Germany only) through MEMMERT service. Requirement: incubator is at its final location and all services have been provided by the purchaser. (Applies only in connection with orders and payments for new units)		K9	K9	K9	

Subject to technical modifications
 □ Standard model, basic specification
 – Not available
 (x) Please specify quantity required after the order number

1) Operating range from 8 °C above ambient
 2) MEMMERT "Celsius" software has been tested for Windows NT 4, 2000, XP and Vista
 3) Option 115V, 50/60 Hz

4) Screws for adjusting height are taken out of the feet and mounted into the subframe
 5) Further temperature distribution surveys at extra cost
 6) Requires Windows 2000 Professional or XP Professional
 7) Service not subject to any discount

Want to know more?
Give Lennox a call on 01 455 2201

OUR PROGRAMME

Universal ovens
Incubators
Hot air sterilisers
Ovens



Vacuum ovens



Peltier-cooled incubators
Compressor-cooled incubators
Cooled incubators



CO₂ incubators



Humidity chambers



Constant climate chambers



Water and oil baths



Your Memmert Partner

Lennox Laboratory Supplies Ltd.
John F. Kennedy Drive,
Naas Road, Dublin 12.
Ireland

email: sales@lennox.ie
web: www.lennox.ie

Upon your request we
are pleased to send you
detailed brochures

memmert[®]

Memmert GmbH + Co. KG
P.O. Box 1720
D-91107 Schwabach
Germany
Tel. +49 (0) 9122 / 925 - 0
Fax +49 (0) 9122 / 145 85
E-Mail: sales@memmert.com

www.memmert.com

Some of the illustrations in this publication include special accessories. We reserve the right of technical modifications. Dimensions are subject to confirmation.